

BW-200 Urine Analyzer

01

BW-200 Urine Analyzer

- [Measurement Principle]** Ultra-high brightness cold light source, reflectometry principle
- [Scope of Application]** 10,11,12 and 14 parameters
- [Test Items]** Glucose (GLU), Bilirubin (BIL), Ketone (KET), Specific gravity (SG), pH, Blood (BLD), Protein (PRO), Urobilinogen (URO), Nitrite (NIT), Leukocyte (LEU), Ascorbic acid (VC), Creatinine (CRE), Calcium (CAL) and Microalbumin (MAL)
- [Test Speed]** 120 samples/hour
- [Record Method]** Thermal print/LCD display
- [Detection Wavelength]** 525nm-660nm
- [Working Environment]** 15°C~35°C, Humidity<85%
- [Power]** 50VA
- [Net Weight]** <3Kg
- [Data Communication]** RS232 interface, can connect with urine sediment analysis system, hospital management system
- [Data Storage]** 2000 test results can be stored, can be queried according to the record number
- [Report Model]** International units or conventional units
- [Printing Method]** Built-in high-speed thermal printer
- [Troubleshooting]** Instrument self-detection with standard calibration strip detection
- [Urine Leakage function]** Automatic urine leakage function, to avoid cross-contamination
- [Volume]** 350MM*285MM*140MM
- [LCD Display]** 240*64 dot matrix, 37MM*130MM



Convenient
Fast
Accurate
Simple



[Characteristic]

- BW-200 is popular with its characters of small and exquisite, convenient and fast and reliable. It is still a classical product until now.

The relentless pursuit of perfect.

