



CE



200mA X-Ray Machine

Brand: Triup Origin: China

Application:
Routine Radiography; Bucky Radiography



TRIUP INTERNATIONAL CORP.

Characteristics:

- Four-sided table-face radiography
- Radiographic clip cassette size: 17"x17 "
- Small focus: 1x1mm², large concentrate: 2x2mm²
- Infinitely varying management of the energy supply voltage V and kilovolt.
- The machine is equipped with subsections and grades mA, kilovolt and S interlocking protective devices.
- The machine is equipped with X-ray filament most protection device, area charging compensation and anode start.
- Power supply voltage, kilovolt photography and mA display separately and are easy to determine.
- No prime and bottom track machine are mounted and easy to operate.
- Transversal upright post movement, long arm movement and-ray tube rotation are all fitted with braking electromagnetism.



Features:

- .Capacity: Two-phase 30KVA.
- .Voltage Frequency: 380V/220V $\pm$ 10% 50Hz $\pm$ 0.5Hz;
- .Internal resistance: 380V: 0.9 Ω , 220V: 0.3 Ω ;
- .Power supply: Adjustable range 380V $\pm$ 10% at 380V; 220 $\pm$ 10% at 220V being adjustable continuously.
- .Time: 0.05-6.3s, 20grades together grading according to R10 coefficient



EMPEX MEDICAL SYSTEM

32, Topkhana Road, Chattogram Samiti Bhaban (1st Floor)
Near Paltan Metro Station, Dhaka-1000 Bangladesh

+8801713594777 +880 2226637337
empextrade@gmail.com empex.com.bd
<https://www.facebook.com/empex.com.bd>



TRIUP INTERNATIONAL CORP.

Address: 8 Floor, FengHuang Culture Building B, No.215,
Jiangdong Middle Rd., Nanjing, China.