



CE



# 300mA X-Ray Machine

Brand: Triup    Origin: China

Application:  
Routine Radiography; Bucky Radiography



TRIUP INTERNATIONAL CORP.

# Characteristics:

- Four-sided table-face radiography
- Radiographic clip cassette size: 17"x17"
- Small focus: 1x1mm<sup>2</sup>, large concentrate: 2x2mm<sup>2</sup>
- Infinitely varying management of the energy supply voltage V and kilovolt.
- The machine is equipped with subsections and grades mA, kilovolt and S interlocking protective devices.
- The machine is equipped with X-ray filament most protection device, area charging compensation and anode start.
- Power supply voltage, kilovolt photography and mA display separately and are easy to determine.
- No prime and bottom track machine are mounted and easy to operate.
- Transversal upright post movement, long arm movement and-ray tube rotation are all fitted with braking electromagnetism.

## Features:

- .Capacity: Two-phase 30KVA.
- .Voltage Frequency: 380V/220V $\pm$ 10% 50Hz $\pm$ 0.5Hz;
- .Internal resistance: 380V: 0.9 $\Omega$ , 220V: 0.3 $\Omega$ ;
- .Power supply: Adjustable range 380V $\pm$ 10% at 380V; 220 $\pm$ 10% at 220V being adjustable continuously.
- .Time: 0.05-6.3s, 20grades together grading according to R10 coefficient



EMPEX MEDICAL SYSTEM



32, Topkhana Road, Chattogram Samiti Bhaban (1st Floor)  
Near Paltan Metro Station, Dhaka-1000 Bangladesh



+8801713594777



+880 2226637337



empextrade@gmail.com



empex.com.bd



<https://www.facebook.com/empex.com.bd>



TRIUP INTERNATIONAL CORP.

Address: 8 Floor, FengHuang Culture Building B, No.215,  
Jiangdong Middle Rd., Nanjing, China.